CHIC AMOUR  
FOLLOW YOUR BLISS

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
**Jojo** is a pair of kegel trainers accompanied by two silky-smooth silicone harnesses. With continued use, Jojo is designed to gradually strengthen your pelvic floor. Each free-floating weight is 28g, with a combined total weight of 56g, and are waterproof for easy cleaning.

## Get to know Jojo

The primary weight can be used individually or combined with the secondary weight using the dual harness.

- 1 Thread the retrieval cord located on the primary weight through either the dual or single harness cord opening.
- 2 Push the primary weight into place, making sure the edges slot firmly into the inner contours of the harness.
- 3 If using the dual harness, push the secondary weight into place, making sure the edges slot firmly into the inner contours of the harness.

Depending on your level of vaginal and pelvic floor fitness, you may wish to start out with the single primary weight and harness and gradually progress to using both weights. The most important thing is that the composition you opt for feels comfortable.


## How to use Jojo

- 1 To facilitate easy insertion, apply a small amount of water-based lubricant on the weight(s) you intend to use. Avoid using too much, as this will make handling the weight(s) difficult.
- 2 While either laying on your back, or standing with one leg raised, gradually introduce the leading end of the weight(s) into your vagina with a firm but gentle push.
- 3 If using the dual harness, insert the secondary weight first. Due to the gap between the two weights, once the secondary is inserted the vaginal muscles may try to push Jojo out before you can insert the primary weight. Take a deep breath and relax.

## Placement

The weight(s) should be positioned so that the lower edge of the primary weight is at least 2 cm inside. The most important thing is that the primary weight rests above your pelvic floor muscle. The exact distance from the vaginal opening to the pelvic floor muscle differs from person to person, but the weight(s) should be positioned comparably to where a tampon would normally be positioned. Always allow the retrieval cord to protrude slightly from the vagina for easier removal.

## Kegel Exercises Tips

For best results, incorporate Kegel exercises into your daily routine and stick to a consistent program.

Depending on your level of pelvic floor fitness, it may take up to 3 months of daily 30-minute Kegel exercising to achieve results.

During the first few weeks of your Kegel routine, you may feel slightly uncomfortable or notice that Jojo tends to slip out after only a few minutes. If this happens, do not worry, this is your body becoming familiar with Jojo. As your PC muscles become stronger, you will become more comfortable and will begin to hold Jojo in place.

When moving about, or even standing in place, the motion and gravity causes the inner free-floating weights to move, creating a subtle dynamic stimulation of the vaginal muscles.

## Caring for Jojo

The life of Jojo will be greatly increased with proper maintenance and care.

- 1 Be sure to wash before first use and after each use.
- 2 Clean using a non-alcohol based antibacterial soap and water, or a dedicated personal massager cleaner.
- 3 Ensure Jojo is completely dry before storing.

## Choosing the proper lubricant

It's important to use the correct lubricant in order to avoid damaging the coating. Jojo is made from non-porous silicone and must be used with a water-based lubricant.

## 90 Day Limited Warranty

If your product fails to work due to any defect in workmanship or materials, we will replace your product with a like or similar item within 90 days of purchase. Simply follow the warranty form procedure located at: **Chicamour.com/warranty**

## Warning

This product is intended for use in kegel training. This product is not a toy. Keep out of reach of children. We are not responsible for the product if the product was used in any other way than the intended use. Store your product on its own and not combined with other devices or personal products. You should refer to the list of materials found on the package to see if the item contains materials that may cause allergic reactions. To avoid injury or aggravation of pre-existing conditions, this device should not be used on swollen, inflamed, infected, irritated areas or skin lacerations. No medical claims are warranted or implied by the use of this product. We do not take any responsibility and we are not liable for any damage caused through use of the product, whether as indirect, special, incidental or consequential damages. Always consult your gynaecologist or midwife for further use. Kegel balls should be avoided during pregnancy and immediately after childbirth, surgery and during menstruation. In case of serious problems such as vaginismus or other similar problems, the use of kegel balls is not recommended.

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